rema<sup>®</sup>  
*TT* 

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Instructions for use

Mode d'emploi

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## Dear Customer,

Thank you for purchasing this product from Dentaaurum. You have chosen a product that will provide you with consistently high quality, reliability and service.

Optimum performance of our dental technology products is based on working precisely to the instructions issued by our Research and Development Department. We have therefore included a clear, detailed summary of each working stage in these instructions for use.

If you have any further questions, our team of dental technicians will be pleased to offer advice.

Simply call: Customer Support Hotline Tel. no. +49 72 31/803-410

Information, instructions for use and safety data sheets for Dentaaurum products are also available on the internet under: [www.dentaaurum.de](http://www.dentaaurum.de).

## Warranty

Dentaaurum guarantees faultless quality of the products manufactured by us. The information in these instructions for use is based upon our own experiences. The user himself is responsible for the correct processing or use of the products. We cannot be held responsible for failures as we have no influence on the processing. However, any claim for compensation only applies to the commercial value of our products.

### Application of the rema® TT test specimen

The test specimen is used to determine the ideal concentration of rema® TT mixing liquid for the production of secondary telescope crowns.

The outer surface of the test specimen has been milled at an optimum angle of 0°. Sometimes laboratory milling machines show deviations from this ideal value (irregular motor, uneven bur). If this is the case for you, please consider the deviation in the actual end result!

#### Tip:

The test specimen can be individualised if desired, i.e.

- by integrating a notch to act as a rotational stop.
- by adding an acrylic or plaster base.
- by re-milling with your own machine to check if your milling machine/bur actually produces a milled angle of 0°.

#### Application

Fabricate a secondary crown in modelling acrylic (e.g. Pattern Resin\*\*), vacuum-formed foil or wax.

Cover the surface with a thin layer of wax.




Invest according to the rema® TT instructions for use.

\* registered trademark of the GC corporation

#### Fit

Too tight → increase the mixing liquid concentration.

Too loose → reduce the mixing liquid concentration.

# Instructions for use – rema® TT

**Phosphate bonded speed investment material for the double crown technique using non-precious metal alloys.**

The special advantages of this investment material are a high degree of material expansion, the freedom of choice concerning pre-heating methods, without having to forfeit the fitting accuracy, and consistently high surface smoothness. In order to benefit from the high quality casting results, it is necessary to observe the following instructions for use carefully.

## Safety Instructions

When using the speed heating method, the furnace door must remain closed for at least 15 minutes after the casting ring has been inserted. (Danger of burning due to spurting flame).

Investment materials contain quartz. Do not inhale the dust! Risk of lung damage (silicosis/lung cancer). Recommendation: Wear protective mask type FFP 2 – EN 149:2001 Cut open the bag using a pair of scissors and avoid generating dust when emptying into the mixing beaker. Fill the empty bag with water before scrunching-up.

Remove dust at the work place using a damp cloth.

In order to avoid dust generation during devestment, soak the ring in water until it is completely moist.

When sandblasting, use an extractor with a fine filter dust bag.

## Delivery

rema® TT Powder	4 kg (25 x 160 g)	REF 106-740-00
rema® TT Mixing Liquid	1000 ml	REF 106-707-00
rema® TT Test Specimen	1 piece	REF 319-740-00

## Service-Life

Powder – 36 months with dry storage

Liquid – 24 months (frost sensitive)

## Material Properties

DIN EN ISO 15912

Type **1**, class **1|2**

Setting begin	6 min*
Compressive strength	9 MPa*
Flowing characteristics	165 mm*
Thermic expansion	1.0 %*

\* Liquid 100 %

# Instructions for Use – rema® TT

Phosphate bonded speed investment material for the double crown technique using non-precious metal alloys.

## Storage

## Wax wetting agent

Only to be used on wax!  
Do not use on modelling acrylic!

## Casting ring system

With the speed heating method,  
only use metal casting rings!

## Mixing ratio

160 g : 33 ml

## Liquid concentration

50 % – 90 %

## Mixing

90 sec. Mixing time under vacuum.

## Working time

4 – 5 min


## Investing

Vibrate gently.

## Conventional heating method

### Additional Information

Observe the ideal working temperature of 20 – 22 °C/8 – 72 °F!  
Only cool the liquid in case of higher room temperatures!  
**Caution!** The total working temperature of powder and liquid is important. For best results, store in a climatic cabinet!

Temperature influences the rate of expansion.

Protect the mixing liquid from frost!  
**Caution with deliveries in winter!**

When using a wax wetting agent, please use Lubrofilm® or Lubrofilm® plus! Please dry the pattern thoroughly!  
Without wax wetting agent, mix thoroughly under vacuum!

Avoid alcohol residues, apply thinly and blow-dry.

For casting ring sizes 3 and 6 please use a dry 2 mm ring liner (REF 127-251-00).

When using silicone or rubber casting rings, these must be soft and elastic.  
Line slightly with Vaseline!

The mixing ratio must be observed accurately!

Use a clean measuring beaker!

### See separate mixing table!

### When using the speed heating method, increase the mixing liquid concentration!

First put liquid, then powder in the mixing beaker.  
Mix by hand thoroughly first.  
Use a clean mixing beaker! No plaster residues!

Mixing speed and mixing blade geometry influence the mixture and therefore the expansion!  
The ideal mixing speed is approx 360 min<sup>-1</sup>.

Measured at room temperature 2 – 22 °C/ 68 – 72 °F.


Cool the liquid at higher temperatures.

Set the vibrator to a low frequency and do not vibrate for too long.



Place the ring liner so that it is flush with the top edge of the ring. Fill the ring completely up to the ring liner. Roughen the upper surface of the investment material after setting!

Setting time: 40 min  
Ideal heating speed: 5 °C/4 °F/min  
Holding time at 250 °C/h  
End temperature: 900 – 950 °C/ 1652 – 1742 °F  
Holding time at end temperature: 40 – 60 min according to ring size


The base of the furnace should be corrugated. Place the ring into the furnace with the opening at the bottom!  
Protect the ring from drying out over the weekend!  
(seal the ring with wax)

### Reduced speed heating

### Speed heating

### Cooling down / devesting Allow to cool down slowly to room temperature.

## Recommended dilution of the mixing liquid

Indication	Premolar	
	Speed heating (metal casting ring)	Conventional (ring-free)
Primary crown/wax	<b>70 %</b>	<b>60 %</b>
Secondary crown (Pattern; Folident)		
Telescope 0° *	<b>70 %</b>	<b>60 %</b>
Conical 2°/4° *	<b>60 %</b>	<b>50 %</b>
For a molar unit, reduce the concentration by approx. 5 %!		
* An optimal fit is achieved, when the secondary crown covers two-thirds of the primary crown after the cast.		

**Tip:** The above values represent the results attained in the Dentaurem test lab.

In individual cases these results can be influenced by different working methods, modelling acrylics and temperatures.

### Ring-free investing !

Setting time: 30 min

Place ring in a pre-heated furnace at 800 °C / 1472 °F!

Residence time: 15 min

Heat up consistently to 900 – 950 °C / 1652 – 1742 °F!

Holding time at end temperature: 60 min

### Only use metal casting rings!

Once the ring has been placed in the furnace, do not open until at least 15 min have elapsed – danger of instant combustion!

Pull the ring liner so that it is flush with the upper edge of the ring. Roughen the surface of the upper investment side with a knife. After 30 min place the ring immediately into the furnace at the end temperature !

(900 – 950 °C / 1652 – 1742 °F)

Holding time: 60 min

Holding time for larger rings (6 – 9): 60 – 90 min

When using plastic parts, always cover with a layer of wax!

Moisten casting rings and carefully deinvest.

Sandblast with glass beads or aluminium oxide at 50 µm.

Do not use a hammer to deinvest!

## Mixing table

Concentration of mixing liquid	160 g : 33 ml	
	Mixing liquid concentration (ml)	Distilled water
90 %	30	3
80 %	26	7
70 %	23	10
60 %	20	13
50 %	17	16

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